5257052104 PQ DE

Resin-based Dental Restorative Material

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CE 1639 REF PE01

## Instructions for Use

Manufacturer

**Tokuyama Dental Corporation**

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

EC Rep.

**Tokuyama Dental Italy S.r.l.**

Via Chizzalunga, 1

36066 Sandrigo, Vicenza, Italy

Tel: +39 0444 659650

Fax: +39 0444 750345

Importer

**Tokuyama Dental Deutschland GmbH**

Am Landwehrbach 5

48341 Altenberge Germany

Tel: +49 2505 938513

Fax: +49 2505 938515

ENGLISH	2	ČESKY	6	DANSK	11	DEUTSCH	16
EESTI KEEL	22	FRANÇAIS	27	ITALIANO	32	LATVISKI	37
LIETUVIŲ K.	42	MAGYAR	47	NEDERLANDS	52	NORSK	57
POLSKI	61	SLOVENČINA	67	SUOMI	71	SVENSKA	76

This product is certified as a medical device in the European Union under the Medical Device Directive 93/42/EEC by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

## ■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE POSTERIOR is a light-cured, radiopaque composite resin for use in posterior restorations. ESTELITE POSTERIOR contains 84% by weight (70% by volume) of silica-zirconia filler (mean particle size: 2  $\mu\text{m}$ , particle size range: 0.1 to 10 $\mu\text{m}$ ). A high filler load offers decreased polymerization shrinkage. The monomer matrix contains Bis-GMA, Triethylene glycol dimethacrylate (TEGDMA) and Bisphenol A polyethoxy methacrylate (Bis-MPEPP).
- 2) ESTELITE POSTERIOR implements Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE POSTERIOR can be shortened by one-third using a light curing unit with a wavelength range of 400 to 500nm (peak: 470nm) compared to our conventional resin composite. Please see the table depicting the relationship between curing time and increment depth (■ INDICATIONS FOR FILLING AND CURING).
- 3) RAP technology facilitates a longer working time for ESTELITE POSTERIOR.
- 4) ESTELITE POSTERIOR is provided in SYRINGE.

## ■ SHADE

PCE, PA1, PA2, PA3.

- Every shade provided for ESTELITE POSTERIOR has adequate translucency designed for posterior restoration.
- PCE serves as a top layer for the multi-shade layering technique.
- Tokuyama Dental Corporation does not provide ESTELITE POSTERIOR Shade Guide. Please make a mockup shade guide by using ESTELITE POSTERIOR.

## ■ INDICATIONS

- Direct posterior restorations including the occlusal surface
- Repair of porcelain/composite

## ■ CONTRAINDICATIONS

ESTELITE POSTERIOR contains methacrylic monomers. DO not use ESTELITE POSTERIOR for patients allergic to or hypersensitive to methacrylic and related monomers.

## ■ PRECAUTIONS

- 1) DO NOT use ESTELITE POSTERIOR for any purpose other than those listed in these instructions. Use ESTELITE POSTERIOR only as directed herein.
- 2) ESTELITE POSTERIOR is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- 3) DO NOT use ESTELITE POSTERIOR if the safety seals are broken or appear to have been tampered with.
- 4) If ESTELITE POSTERIOR causes an allergic reaction or oversensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE POSTERIOR to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE POSTERIOR comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of ESTELITE POSTERIOR with eyes, mucosal membrane, skin and clothing.
  - If ESTELITE POSTERIOR comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
  - If ESTELITE POSTERIOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.
  - If ESTELITE POSTERIOR comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
  - Instruct the patient to rinse his mouth immediately after treatment.
- 7) ESTELITE POSTERIOR should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of ESTELITE POSTERIOR, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

## ■ INDICATIONS FOR FILLING AND CURING

Be sure to light-cure ESTELITE POSTERIOR extraorally and check the time needed for complete hardening of ESTELITE POSTERIOR with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

## Relationship between curing time and increment depth:

Light Type	Intensity (mW/cm <sup>2</sup> )	Curing Time (seconds)	Increment Depth (mm) <sup>1)</sup>
Halogen	400	10	1.9
		20	2.3
	600	10	2.0
		20	2.4
	800	5	1.8
		10	2.4
LED	260	10	1.9
		20	2.0
	900	5	1.8
		10	2.0
Plasma Arc	950	3	2.1

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 “depth of cure” of ISO4049.

## ■ STORAGE

- 1) Store ESTELITE POSTERIOR at a temperature between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use ESTELITE POSTERIOR after the indicated date of expiration on the syringe or package.

## ■ DISPOSAL

To safely dispose of excess ESTELITE POSTERIOR, extrude unused portion from SYRINGE and light-cure before disposal.

## ■ CLINICAL PROCEDURE

### 1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

### 2. Shade Selection

Select the appropriate shade of ESTELITE POSTERIOR prior to placing the rubber dam.

- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist

teeth.

- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

### **3. Isolation**

A rubber dam is the preferred method to isolate the tooth from contamination.

### **4. Cavity Preparation**

Prepare the cavity and rinse with water. Add chamfers to the margins of posterior preparations (class I, II) because chamfers assist in diminishing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

### **5. Pulp Protection**

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL BASED MATERIALS** to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE POSTERIOR.

### **6. Bonding system**

Apply a bonding system according to manufacturer's instructions.

### **7. Dispensing**

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. Re-cap the syringe immediately after dispensing.

- **DO NOT** apply unnecessary force to the syringe immediately after removing from the refrigerator.

### **8. Filling and contouring**

Fill the cavity incrementally using a suitable placing instrument. Increments should not exceed the indicated depth (please refer to the aforementioned table).

- **DO NOT** mix ESTELITE POSTERIOR with other brands of resin composite or other shades of ESTELITE POSTERIOR to avoid incomplete cure or entrapment of air bubbles.

### **9. Curing**

Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

## 10. Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To diminish the demarcation between the composite and the natural tooth, rotate a 12 fluted carbide bur slowly along the margins without water. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGo™ (DENTSPLY/Caulk), Sof-Lex™ (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINE™ Hybrid Diamond (COSMEDENT, Inc.).

\* PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

- Instruct the patient adequately on the maintenance of the restoration after treatment.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**IMPORTANT NOTE:** The manufacturer of ESTELITE POSTERIOR is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use. Product specifications of ESTELITE POSTERIOR are subjected to change without notification. When the product specification change, the instruction and precautions may change also.

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\*PoGo, Sof-Lex, Identoflex och D-FINE är varumärken som inte tillhör Tokuyama Dental Corporation.