Light-Cured Dental Adhesive

TOKUYAMA EE-BONDCE₁₆₃₉ REF AD02**Instructions for Use**

Manufacturer

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

EC Rep.

Tokuyama Dental Italy S.r.l.

Via Chizzalunga, 1

36066 Sandrigo, Vicenza, Italy

Tel: +39 0444 659650

Fax: +39 0444 750345

Importer

Tokuyama Dental Deutschland GmbH

Am Landwehrbach 5

48341 Altenberge Germany

Tel: +49 2505 938513

Fax: +49 2505 938515

ENGLISH	2	ČESKY	8	DANSK	14	DEUTSCH	20
EESTI KEEL	26	FRANÇAIS	33	ITALIANO	39	LATVISKI	46
LIETUVIŲ K.	52	MAGYAR	58	NEDERLANDS	64	NORSK	70
POLSKI	76	SLOVENČINA	83	SUOMI	89	SVENSKA	95

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1. TOKUYAMA EE-BOND is a light-cured fluoride releasing dental adhesive system which contains TOKUYAMA ETCHING GEL HV (etching agent for enamel etching technique) and EE-BOND (one component light-cured bonding agent for dentin and enamel).
2. Enamel etching technique allows EE-BOND to form a durable bonding layer on dentin and enamel due to good penetration of EE-BOND into both tooth structures. TOKUYAMA EE-BOND exhibits excellent adhesive properties and marginal integrity to enamel and dentin when used in combination with light- or dual-cured composite materials. Its excellent marginal integrity to uncut enamel will enhance esthetic restorations.
3. A light-curing unit having the camphorquinone (CQ) wavelength range (peak: 470nm, spectrum: 400 to 500nm) can be used for curing EE-BOND.
4. TOKUYAMA ETCHING GEL HV contains 39 wt% Phosphoric acid, Purified water, Thickener, Colorant. EE-BOND contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate, 2-Hydroxyethyl methacrylate (HEMA), Camphorquinone, and solvent.
5. TOKUYAMA ETCHING GEL HV and EE-BOND are dispensed in syringe and bottle, respectively.

■ INDICATIONS

Bonding of light- or dual-cured composite material to:

- cut/uncut enamel,
- cut/uncut dentin,
- fractured porcelain/composite repair.

■ CONTRAINDICATIONS

1. DO NOT use TOKUYAMA EE-BOND for patients allergic to or hypersensitive to acids, methacrylic monomers, related monomers, or organic solvents.
2. DO NOT mix TOKUYAMA ETCHING GEL HV and sodium hypochlorite, mixing will generate Chloride gas which is harmful.

■ PRECAUTIONS

- 1) DO NOT use TOKUYAMA EE-BOND for any purpose other than those listed in these instructions. Use TOKUYAMA EE-BOND only as directed herein.
- 2) TOKUYAMA EE-BOND is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- 3) DO NOT use TOKUYAMA EE-BOND if the safety seals are broken or appear to have been tampered with.
- 4) If TOKUYAMA EE-BOND causes an allergic reaction or oversensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling TOKUYAMA EE-BOND to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If TOKUYAMA ETCHING GEL HV or EE-BOND come in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of TOKUYAMA ETCHING GEL HV or EE-BOND with eyes, mucosal membrane, skin and clothing.
 - If TOKUYAMA ETCHING GEL HV or EE-BOND comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If TOKUYAMA ETCHING GEL HV comes in contact with the mucosal membrane and skin, wipe the affected area immediately, and thoroughly flush with water and immediately contact a physician.
 - If EE-BOND comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration. Affected areas may whiten from protein coagulation, but such whitening should disappear within 24 hours. If such whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If EE-BOND comes in contact with the skin, immediately saturate the area with alcohol-soaked cotton swab or gauze and flush with water.
 - If TOKUYAMA ETCHING GEL HV or EE-BOND comes in contact with clothing immediately saturate the area with alcohol-soaked cotton swab or gauze and flush with water.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) TOKUYAMA EE-BOND should not be ingested or aspirated. Ingestion or aspiration may

cause serious injury.

- 8) To avoid the unintentional ingestion of TOKUYAMA EE-BOND, do not leave it unsupervised within the reach of patients and children.
- 9) DO NOT expose EE-BOND or its vapor to open flame because EE-BOND is flammable.
- 10) To avoid cross contamination, DO NOT reuse the disposable tip and disposable applicator included in the TOKUYAMA EE-BOND package.
- 11) Clean the dispensing well included in the TOKUYAMA EE-BOND package thoroughly with alcohol after each use.
- 12) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.

■ PRECAUTIONS FOR MEDICAMENTS AND MATERIALS

- 1) Some materials and medicaments (hemostatic material) inhibit adhesion of TOKUYAMA EE-BOND for an extended period, even after meticulous cleansing with water. DO NOT USE products which contain:
 - eugenol,
 - hydrogen peroxide,
 - sodium hypochlorite,
 - diammine silver fluoride [molecular formula: $\text{Ag}(\text{NH}_3)_2\text{F}$],
 - phenols such as parachlorophenol, guaiacol, phenol,
 - aluminum chloride,
 - ferric sulfate,
 - aluminum sulfate,
 - epinephrine.
- 2) To avoid the mixing of TOKUYAMA ETCHING GEL HV and sodium hypochlorite (sodium hypochlorite included materials and medicaments) when using both in the same procedure, thoroughly rinse before using either.

■ STORAGE

- 1) Store TOKUYAMA ETCHING GEL HV at temperatures between 0 to 25°C (32 to 77°F)
- 2) Store EE-BOND under refrigeration at temperatures between 0 to 10 °C (32 to 50°F).
- 3) Keep away from heat, direct sunlight, sparks and open flames.
- 4) DO NOT use after the indicated expiration date on the syringe/bottle/package.

■ DISPOSAL

Unused TOKUYAMA ETCHING GEL HV and EE-BOND are to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■ CLINICAL PROCEDURES

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation

The rubber dam is the preferred method of isolation.

3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply TOKUYAMA ETCHING GEL HV for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling reagent according to its manufacturer's instructions.

4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- The substances listed below, which inhibit curing of EE-BOND, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of TOKUYAMA ETCHING GEL HV for 2 to 3 seconds before application:

- 1) Oil mist from handpiece,
- 2) Saliva, blood and exudates.

5. Pulp protection

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of EE-BOND.

6. Enamel etching

Attach the Disposable Tip after removing the cap of TOKUYAMA ETCHING GEL HV. Verify the flow of TOKUYAMA ETCHING GEL HV prior to applying intra-orally. Apply

TOKUYAMA ETCHING GEL HV (39 wt% phosphoric acid etching gel) only to the uncut enamel surrounding margin of prepared cavity and leave TOKUYAMA ETCHING GEL HV in place for 5 seconds. Rinse the etched surface thoroughly (at least 5 seconds) with water, and then dry with mild air. Remove the Disposable Tip from TOKUYAMA ETCHING GEL HV and replace the syringe cap.

- Overfilling of composite materials onto uncut enamel not being etched could cause marginal microleakage and discoloration.
- Cut enamel does not need to be etched. Applying TOKUYAMA ETCHING GEL HV to cut enamel will not improve or impair the adhesive properties of EE-BOND toward cut enamel.
- Applying TOKUYAMA ETCHING GEL HV to dentin may reduce the bond strength of EE-BOND to dentin.

7. Dispensing of EE-BOND

Open the bottle cap of EE-BOND and dispense one or two drops of EE-BOND onto the Dispensing Well. Close the bottle cap tightly immediately after dispensing.

- Wipe off the excess adhesive on the tip of nozzle before closing.
- Do not mix the adhesive with other brands of primers or adhesives.

8. Application of EE-BOND

Apply EE-BOND to the cavity walls, margin and surrounding etched uncut enamel using the disposable applicator. Be sure not to leave out any areas where the EE-BOND should be applied. Leave undisturbed for 10 seconds after the end of application.

- Protect the dispensed EE-BOND and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the EE-BOND application within 5 minutes after dispensing due to the volatile alcohol content.
- In case of multiple restorations, ensure individual application time for each restoration.
- If saliva, blood, or other fluids contaminate the applied EE-BOND, thoroughly rinse the cavity with water, dry and re-apply EE-BOND.
- Do not rinse the applied EE-BOND with water except in unintentional contamination.

9. Air Dry

Using an oil-free air/water syringe, apply weak air flow to the adhesive surface, continuing until the runny EE-BOND stays in the same position without any motion (usually for 5 seconds). Finish by using a mild air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the EE-BOND.

- If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.

- To prevent spattering, refer to the following tips:
 - 1) To avoid inadvertent strong air;
 - a) Be sure to begin weak air flow outside the mouth,
 - b) Direct the air flow at the EE-BOND surface.
 - 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror refraction technique is also useful in reducing the air flow.
 - 3) If the EE-BOND pools on the cavity floor or the cavosurface angle and is too thick to air thin, blot the excess with a new disposable applicator before weak air application

10. Light-cure

Light-cure the surface for 10 seconds or more, keeping the curing light tip within a distance of 2 mm from the surface. If the cavity is too large or too distant (e.g. MOD), divide the area into segments and light-cure each segment individually.

- Confirm that the light-curing unit has sufficient intensity ($>300\text{mW}/\text{cm}^2$) before using. Note that using a cracked light guide will lower the intensity.

11. Light-cured Composite

Restore with light-cured composite resin according to its manufacturer's instructions. The overfilled composite resin should be thoroughly finished and polished.

- When dual-cured composite resins are placed into a cavity, the first increment must be light-cured using a layering technique.
- Do not use self-cured composite resins here because the Phosphoric acid monomer contained in EE-BOND may interfere with the curing of the self-cured resins, leading to their premature detachment.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The manufacturer of the TOKUYAMA EE-BOND is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Product specifications of the TOKUYAMA EE-BOND are subject to change without notification. When the product specifications change, the instructions and precautions may change also.

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