Resin-based Dental Restorative Material

ESTELITE ASTERIA SYRINGEESTELITE ASTERIA PLT

Instructions for Use

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

| | | | | | | | |
|--------------------|----|-------------------|----|-------------------|----|-----------------|-----|
| ENGLISH | 2 | ČESKÝ | 8 | DANSK | 14 | DEUTSCH | 21 |
| EESTI KEEL | 28 | FRANÇAIS | 35 | ITALIANO | 42 | LATVISKI | 49 |
| LIETUVIŲ K. | 55 | MAGYAR | 62 | NEDERLANDS | 69 | NORSK | 75 |
| POLSKI | 81 | SLOVENČINA | 89 | SUOMI | 95 | SVENSKA | 101 |

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE ASTERIA is a light-cured, radiopaque composite resin for use in anterior and posterior restorations. It is indicated for all carious lesion classes, including minimal or no-preparation type restorations. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. ESTELITE ASTERIA contains 82% by weight (71% by volume) silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. Every inorganic filler contained in ESTELITE ASTERIA is a spherical filler (mean particle size: 200 nm, particle size range: 100 to 300 nm) that enables excellent gloss retention and wear resistance. ESTELITE ASTERIA contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), 1,6-bis(methacryl-ethylxycarbonylamino)trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinol, Dibutyl hydroxyl toluene, and UV absorber.
- 2) ESTELITE ASTERIA incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time (90 sec.). Please see the table depicting the relationship between curing time and increment depth (■ INDICATIONS FOR FILLING AND CURING).
- 3) ESTELITE ASTERIA is provided either in SYRINGE or Pre Loaded Tip(PLT).

■ SHADE

ESTELITE ASTERIA contains 7 Body shades and 5 Enamel shades. The ESTELITE ASTERIA shading system is designed to mimic a wide shade range of the natural dentition by layering the appropriate shades selected from 2 categories: Body and Enamel.

| Component | Shade | Tip |
|-----------|------------------------------------|---|
| Body | A1B, A2B, A3B, A3.5B, A4B, B3B, BL | Body shades are designed for reconstructing the dentin layer. The Body shades <u>should cover all enamel margins except the incisal area.</u> |
| Enamel | Natural Enamel (NE) | NE is recommended to restore translucency in the incisal area in most cases. |
| | White Enamel (WE) | WE is recommended for the proximal wall. WE is <u>suggested as an alternative to NE in whiter cases.</u> |
| | Yellow Enamel (YE) | YE is designed to mimic discolored enamel. |
| | Trans Enamel (TE) | TE is the most translucent in ESTELITE ASTERIA. This shade is suggested as an alternative to NE in <u>highly translucent cases.</u> |
| | Occlusal Enamel (OcE) | OcE is recommended for the occlusal surface. OcE has exceptional sculptability to shape occlusal cusps and fissures. |

■ INDICATIONS

- Direct anterior and posterior restorations including occlusal surfaces
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

■ CONTRAINDICATIONS

ESTELITE ASTERIA contains methacrylic monomers. DO NOT use ESTELITE ASTERIA for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

■ PRECAUTIONS

- 1) DO NOT use ESTELITE ASTERIA for any purposes other than those listed in these instructions. Use ESTELITE ASTERIA only as directed herein.
- 2) ESTELITE ASTERIA is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use ESTELITE ASTERIA if the safety seals are broken or appear to have been tampered with.

- 4) If ESTELITE ASTERIA causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTELITE ASTERIA to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE ASTERIA comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of ESTELITE ASTERIA with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE ASTERIA comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE ASTERIA comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE ASTERIA comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) ESTELITE ASTERIA should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of ESTELITE ASTERIA, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments and brushes with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

INDICATIONS FOR FILLING AND CURING

ESTELITE ASTERIA is designed to be cured by either a halogen, LED or Plasma Arc light with a wavelength of 400-500 nm. Be sure to light-cure ESTELITE ASTERIA extra-orally and check the time needed for complete hardening of ESTELITE ASTERIA with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

| Light type | Intensity (mW/cm ²) | Curing time (Seconds) | Increment depth (mm) ¹⁾ | |
|------------|---------------------------------|-----------------------|---|-----------------|
| | | | A1B, A2B, A3B, A3.5B, B3B, BL NE, WE, YE, TE, OcE | A4B |
| Halogen | 400 | 10 | 1.7 | – ²⁾ |
| | | 20 | 2.1 | 1.8 |
| | 600 | 10 | 2.0 | 1.7 |
| | | 20 | 2.4 | 2.0 |
| | 800 | 10 | 2.1 | 1.8 |
| | | 20 | 2.5 | 2.1 |
| LED | 400 | 10 | 1.7 | – ²⁾ |
| | | 20 | 2.0 | 1.7 |
| | 600 | 10 | 1.8 | 1.6 |
| | | 20 | 2.2 | 1.9 |
| | 1000 | 10 | 2.0 | 1.8 |
| | | | | |
| Plasma Arc | 950 | 6 | 1.9 | 1.8 |

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 “depth of cure” of ISO4049.

2) Curing time must be at least 20 seconds.

■ SPECIAL NOTES FOR THE USE OF PLT

- 1) PLTs are designed for single patient use only. To avoid cross infection, Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
- 2) Dispensers are not provided for ESTELITE ASTERIA. Use a dispenser that fits the PLT of ESTELITE ASTERIA. For proper use and proper disinfection, see manufacturer’s instructions.
- 3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■ STORAGE

- 1) Store ESTELITE ASTERIA at temperatures between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use ESTELITE ASTERIA after the indicated date of expiration on the syringe or PLT package.

■ DISPOSAL

To safely dispose of excess ESTELITE ASTERIA, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

Follow local instructions for disposal.

Refer to the packaging waste disposal guide shown on the end of this IFU.

■ CUSTOM SHADE GUIDE

See the instructions for ESTELITE ASTERIA Custom Shade Guide prior use.

■ CLINICAL PROCEDURE

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste, and then rinse with water.

2. Shade Selection

Select the appropriate shade of ESTELITE ASTERIA using the custom shade guide prepared in advance.

- Complete the shade selection within 5 minutes; teeth become whiter when dehydrated.
- Lightness (color value) is the primary importance on shade selection.
- In the case of whitened teeth, select the shade a few weeks later after the completion of whitening. Whitened teeth tend to become slightly darker with the elapse of time.
- Please see the table depicting the shade and tip (■ SHADE).

3. Isolation

A rubber dam is the preferred method of isolation.

4. Cavity Preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), bevels assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.

- Scalloped (undulant) bevel may be desirable depending on the esthetics and retention.
- In the case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to

manufacturer's instructions.

5. Pulp Protection

Protect the pulp appropriately, if the cavity is in close proximity to the pulp. In the case of pulp exposure, apply calcium hydroxide for pulp capping. **DO NOT USE EUGENOL-BASED MATERIALS** to protect the pulp as these materials will inhibit curing of ESTELITE ASTERIA.

6. Bonding System

Apply a bonding system according to its manufacturer's instructions.

7. Dispensing

7-1. PLT

Please read ■ **SPECIAL NOTES FOR THE USE OF PLT** before use.

- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into cavity or indirectly on the mixing pad.

7-2. SYRINGE

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.

- **DO NOT** apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring

Fill the cavity incrementally. Increments should not exceed the indicated curing depth (please refer to the aforementioned table).

- In case that characterization is required, tints (such as ESTELITE COLOR) can be used. Apply the tint according to manufacturer's instructions.
- **DO NOT** mix ESTELITE ASTERIA with other brands of composite resin to avoid incomplete cure or entrapment of air bubbles. In order to avoid air bubbles entrapment, **DO NOT** mix ESTELITE ASTERIA with another shade of the paste.

9. Curing

Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow its manufacturer's instructions.

10. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing

diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

- Instruct the patient adequately on the maintenance of the restoration after treatment.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

ENGLISH

Packaging waste disposal guide

| | |
|---------------------------------|--|
| SYRINGE (cap and body) | Dispose of as industrial waste. |
| PLT (cap and body) | Dispose of as industrial waste. |
| Packaging bag for PLT | PP5 (Polypropylene): Recycling. Follow local instructions for disposal. |
| Case for PLT | PP5 (Polypropylene): Recycling. Follow local instructions for disposal. |
| Case for SYRINGE (cap and body) | SBR7 (Styrene-butadiene rubber): Recycling. Follow local instructions for disposal. |
| Case partation for SYRINGE | PET1 (Polyethylene terephthalate): Recycling. Follow local instructions for disposal. |
| Paperbox and/or partation | PAP21 (Non-corrugated fibreboard): Recycling. Follow local instructions for disposal. |
| IFU and/or guide | PAP22 (Paper): Recycling. Follow local instructions for disposal. |

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