

# SAFETY DATA SHEET

## ESTELITE Σ QUICK

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product name ESTELITE Σ QUICK

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Identified uses Resin-based Dental Restorative Material. For dental professionals only.

#### 1.3. Details of the supplier of the safety data sheet

##### Supplier

Tokuyama Dental Italy S.r.l.  
Via Chizzalunga, 1,  
36066 Sandrigo, Vicenza, Italy  
TEL: +39-0444-659650  
FAX: +39-0444-750345

##### Contact person

<http://www.tokuyama-dental.com/tdc/contact.html>

##### Manufacturer

Tokuyama Dental Corporation  
38-9, Taitou 1-chome, Taitou-ku, Tokyo  
110-0016, Japan  
TEL: +81-3-3835-2261  
FAX: +81-3-3835-2265

#### 1.4. Emergency telephone number

Emergency telephone Regional Medicines and Poisons Information Centre NI,  
Pharmacy Department  
Royal Hospital Suite, Grosvenor Road, Belfast  
Emergency telephone: 844 892 0111  
Telephone: +44 28 90 63 2032  
Fax: +44 28 90 24 8030  
E-mail address: [nirdic.nirdic@belfasttrust.hscni.net](mailto:nirdic.nirdic@belfasttrust.hscni.net)

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### Classification (EC 1272/2008)

Physical hazards Not Classified  
Health hazards Skin Sens. 1 - H317  
Environmental hazards Not Classified

#### 2.2. Label elements

##### Pictogram



##### Signal word

Warning

##### Hazard statements

H317 May cause an allergic skin reaction.

## ESTELITE Σ QUICK

<b>Precautionary statements</b>	<p>P261 Avoid breathing vapours.</p> <p>P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.</p> <p>P302+P352 IF ON SKIN: Wash with plenty of water.</p> <p>P501 Dispose of contents/ container in accordance with local regulations.</p>
<b>Contains</b>	(1-METHYLETHYLIDENE)BIS[4,1-PHENYLENEOXY(2-HYDROXY-3,1-PROPANEDIYL)] BISMETHACRYLATE, 2,2'-ETHYLENEDIOXYDIETHYL DIMETHACRYLATE, MEQUINOL
<b>Supplementary precautionary statements</b>	<p>P272 Contaminated work clothing should not be allowed out of the workplace.</p> <p>P321 Specific treatment (see medical advice on this label).</p> <p>P333+P313 If skin irritation or rash occurs: Get medical advice/ attention.</p> <p>P362+P364 Take off contaminated clothing and wash it before reuse.</p>

### 2.3. Other hazards

This product does not contain any substances classified as PBT or vPvB.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

<b>(1-METHYLETHYLIDENE)BIS[4,1-PHENYLENEOXY(2-HYDROXY-3,1-PROPANEDIYL)] BISMETHACRYLATE</b> <span style="float: right;"><b>10-30%</b></span>
CAS number: 1565-94-2                      EC number: 216-367-7
<b>Classification</b> Skin Sens. 1 - H317 Aquatic Chronic 3 - H412
<b>2,2'-ETHYLENEDIOXYDIETHYL DIMETHACRYLATE</b> <span style="float: right;"><b>5-10%</b></span>
CAS number: 109-16-0                      EC number: 203-652-6
<b>Classification</b> Skin Irrit. 2 - H315 Eye Irrit. 2 - H319 Skin Sens. 1 - H317
<b>TITANIUM DIOXIDE</b> <span style="float: right;"><b>&lt;1%</b></span>
CAS number: 13463-67-7                      EC number: 236-675-5
<b>Classification</b> Acute Tox. 4 - H332 Carc. 2 - H351
<b>MEQUINOL</b> <span style="float: right;"><b>&lt;1%</b></span>
CAS number: 150-76-5                      EC number: 205-769-8
<b>Classification</b> Acute Tox. 4 - H302 Eye Irrit. 2 - H319 Skin Sens. 1 - H317

# ESTELITE Σ QUICK

<b>2,6-DI-tert-BUTYL-p-CRESOL</b>	<b>&lt;1%</b>
CAS number: 128-37-0	EC number: 204-881-4
M factor (Chronic) = 1	
<b>Classification</b> Aquatic Chronic 1 - H410	

The full text for all hazard statements is displayed in Section 16.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

<b>General information</b>	Move affected person to fresh air at once. Get medical attention if any discomfort continues.
<b>Inhalation</b>	Move affected person to fresh air at once.
<b>Ingestion</b>	Try to induce vomiting. Get medical attention if any discomfort continues.
<b>Skin contact</b>	Wash skin thoroughly with soap and water.
<b>Eye contact</b>	Make sure to remove any contact lenses from the eyes before rinsing. Promptly wash eyes with plenty of water while lifting the eye lids. Continue to rinse for at least 15 minutes. Get medical attention if irritation persists after washing.

### 4.2. Most important symptoms and effects, both acute and delayed

<b>Skin contact</b>	Allergic rash.
---------------------	----------------

### 4.3. Indication of any immediate medical attention and special treatment needed

<b>Notes for the doctor</b>	No specific recommendations.
-----------------------------	------------------------------

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

<b>Suitable extinguishing media</b>	Extinguish with the following media: Foam. Carbon dioxide or dry powder.
-------------------------------------	--

### 5.2. Special hazards arising from the substance or mixture

<b>Specific hazards</b>	When heated and in case of fire, irritating vapours/gases may be formed.
<b>Hazardous combustion products</b>	No known hazardous decomposition products.

### 5.3. Advice for firefighters

<b>Protective actions during firefighting</b>	Avoid breathing fire gases or vapours.
<b>Special protective equipment for firefighters</b>	Wear chemical protective suit.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

<b>Personal precautions</b>	Wear protective clothing as described in Section 8 of this safety data sheet.
-----------------------------	---

### 6.2. Environmental precautions

<b>Environmental precautions</b>	No specific recommendations.
----------------------------------	------------------------------

### 6.3. Methods and material for containment and cleaning up

# ESTELITE Σ QUICK

**Methods for cleaning up** Eliminate all sources of ignition. No smoking, sparks, flames or other sources of ignition near spillage. Provide adequate ventilation. Absorb in vermiculite, dry sand or earth and place into containers.

## 6.4. Reference to other sections

**Reference to other sections** For waste disposal, see Section 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

**Usage precautions** All handling should only take place in well-ventilated areas.

### 7.2. Conditions for safe storage, including any incompatibilities

**Storage precautions** Store in a cool and well-ventilated place. Store at temperatures between 0°C/25°F and 32°C/77°F. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid exposure to high temperatures or direct sunlight.

**Storage class** Chemical storage.

### 7.3. Specific end use(s)

**Specific end use(s)** The identified uses for this product are detailed in Section 1.2.

## SECTION 8: Exposure Controls/personal protection

### 8.1. Control parameters

#### Occupational exposure limits

##### TITANIUM DIOXIDE

Long-term exposure limit (8-hour TWA): WEL 4 mg/m<sup>3</sup> respirable dust

Long-term exposure limit (8-hour TWA): WEL 10 mg/m<sup>3</sup> inhalable dust

##### 2,6-DI-tert-BUTYL-p-CRESOL

Long-term exposure limit (8-hour TWA): WEL 10 mg/m<sup>3</sup>

WEL = Workplace Exposure Limit

### 8.2. Exposure controls

**Appropriate engineering controls** Provide adequate general and local exhaust ventilation.

**Eye/face protection** The following protection should be worn: Chemical splash goggles.

**Hand protection** Wear protective gloves.

**Other skin and body protection** Wear appropriate clothing to prevent any possibility of skin contact.

**Hygiene measures** Do not smoke in work area. Wash at the end of each work shift and before eating, smoking and using the toilet. Wash promptly with soap and water if skin becomes contaminated. Promptly remove any clothing that becomes contaminated. Use appropriate skin cream to prevent drying of skin. When using do not eat, drink or smoke.

## SECTION 9: Physical and Chemical Properties

### 9.1. Information on basic physical and chemical properties

**Appearance** Paste

**Colour** Yellow-white

**Odour** Not available.

## ESTELITE Σ QUICK

Odour threshold	Not available.
pH	pH (concentrated solution): 7
Melting point	Not available.
Initial boiling point and range	Not available.
Flash point	Not available.
Evaporation rate	Not available.
Evaporation factor	Not available.
Flammability (solid, gas)	Not available.
Upper/lower flammability or explosive limits	Not applicable.
Other flammability	Not available.
Vapour pressure	Not available.
Vapour density	Not available.
Relative density	1.9
Solubility(ies)	Not available.
Partition coefficient	Not available.
Auto-ignition temperature	Not available.
Decomposition Temperature	Not available.
Viscosity	Not available.
Explosive properties	Not explosive.
Oxidising properties	Not available.

### 9.2. Other information

Other information Not available.

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

Reactivity May polymerise. See the other subsections of this section for further details.

### 10.2. Chemical stability

Stability Stable at normal ambient temperatures and when used as recommended.

### 10.3. Possibility of hazardous reactions

Possibility of hazardous reactions Not known.

### 10.4. Conditions to avoid

Conditions to avoid Keep away from heat, sparks and open flame. Protect against direct sunlight.

### 10.5. Incompatible materials

Materials to avoid Strong oxidising agents.

### 10.6. Hazardous decomposition products

## ESTELITE Σ QUICK

**Hazardous decomposition products** Heating may generate the following products: Toxic gases or vapours.

### SECTION 11: Toxicological information

#### 11.1. Information on toxicological effects

**Routes of exposure** See section 4.2. Skin , Eyes , Ingestion , Inhalation

#### Acute toxicity - oral

**Notes (oral LD<sub>50</sub>)** Not available.

#### Acute toxicity - dermal

**Notes (dermal LD<sub>50</sub>)** Not available.

#### Acute toxicity - inhalation

**Notes (inhalation LC<sub>50</sub>)** Not available.

#### Skin corrosion/irritation

**Skin corrosion/irritation** Not available.

#### **Animal data**

Not available.

#### Serious eye damage/irritation

**Serious eye damage/irritation** Not available.

#### Respiratory sensitisation

**Respiratory sensitisation** Not available.

#### Skin sensitisation

**Skin sensitisation** Sensitising.

#### Germ cell mutagenicity

**Genotoxicity - in vitro** Not available.

**Genotoxicity - in vivo** Not available.

#### Carcinogenicity

**Carcinogenicity** Not available.

**IARC carcinogenicity** Some of the ingredients are listed or exempt.

**NTP carcinogenicity** Some of the ingredients are listed or exempt.

#### Reproductive toxicity

**Reproductive toxicity - fertility** Not available.

**Reproductive toxicity - development** Not available.

#### Specific target organ toxicity - single exposure

**STOT - single exposure** Not available.

#### Specific target organ toxicity - repeated exposure

**STOT - repeated exposure** Not available.

#### Aspiration hazard

**Aspiration hazard** Not available.

### TITANIUM DIOXIDE

## ESTELITE Σ QUICK

### Acute toxicity - inhalation

ATE inhalation 1.5  
(dusts/mists mg/l)

### Carcinogenicity

IARC carcinogenicity IARC Group 2B Possibly carcinogenic to humans.

NTP carcinogenicity Reasonably anticipated to be a human carcinogen.

OSHA Carcinogenicity Not listed.

### 2,6-DI-tert-BUTYL-p-CRESOL

### Carcinogenicity

IARC carcinogenicity IARC Group 3 Not classifiable as to its carcinogenicity to humans.

NTP carcinogenicity Reasonably anticipated to be a human carcinogen.

OSHA Carcinogenicity Not listed.

## SECTION 12: Ecological Information

### 12.1. Toxicity

Acute toxicity - fish Not available.

Acute toxicity - aquatic invertebrates Not available.

Acute toxicity - aquatic plants Not available.

Chronic toxicity - fish early life stage Not available.

Short term toxicity - embryo and sac fry stages Not available.

Chronic toxicity - aquatic invertebrates Not available.

### 12.2. Persistence and degradability

Persistence and degradability No data available.

### 12.3. Bioaccumulative potential

Bioaccumulative potential No data available on bioaccumulation.

Partition coefficient Not available.

### 12.4. Mobility in soil

Mobility No information available.

### 12.5. Results of PBT and vPvB assessment

Results of PBT and vPvB assessment This product does not contain any substances classified as PBT or vPvB.

### 12.6. Other adverse effects

Other adverse effects Not known.

## SECTION 13: Disposal considerations

# ESTELITE Σ QUICK

## 13.1. Waste treatment methods

**Disposal methods** Dispose of waste to licensed waste disposal site in accordance with the requirements of the local Waste Disposal Authority.

## SECTION 14: Transport information

**General** The product is not covered by international regulations on the transport of dangerous goods (IMDG, IATA, ADR/RID).

### 14.1. UN number

Not applicable.

### 14.2. UN proper shipping name

Not applicable.

### 14.3. Transport hazard class(es)

Not applicable.

### 14.4. Packing group

Not applicable.

### 14.5. Environmental hazards

**Environmentally hazardous substance/marine pollutant**

No.

### 14.6. Special precautions for user

Not applicable.

### 14.7. Transport in bulk according to Annex II of MARPOL and the IBC Code

**Transport in bulk according to** Not applicable.

**Annex II of MARPOL 73/78  
and the IBC Code**

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

**EU legislation** (EC) No 1907/2006 (REACH).  
(EC) No 1272/2008 (CLP).  
(EU) No 2015/830.

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out.

## SECTION 16: Other information

**Classification procedures according to Regulation (EC) 1272/2008** Skin Sens. 1 - H317: Calculation method.

**Training advice** Ensure operatives are trained to minimise exposure. Only trained personnel should use this material.

**Revision comments**

**Revision date** 15/09/2018



## ESTELITE Σ QUICK

<b>Revision</b>	3
<b>Hazard statements in full</b>	H302 Harmful if swallowed. H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H332 Harmful if inhaled. H351 Suspected of causing cancer. H410 Very toxic to aquatic life with long lasting effects. H412 Harmful to aquatic life with long lasting effects.

The information which is contained in this document is based on available data. However, as such has been obtained from various sources, including independent laboratories, it is given without warranty or representation that is complete, accurate and can be relied upon. Tokuyama Dental Corp. has not attempted to conceal in any way the deleterious aspects of the product listed herein, but makes no warranty as to such information.